
Samsung Medison is a manufacturer of medical devices. Founded in 1985, the company sells cutting-edge medical devices, including diagnostic ultrasound, digital X-ray, and blood analysis machines, and is known around the world for its R&D capabilities and advanced technologies. In 2011, it became an affiliate of Samsung Electronics.

- * This product, along with its various features, options and transducers, is not currently available in all countries. Due to regulatory reasons, its future availability cannot be guaranteed. Please contact your local sales network for further details.
- * Beyond Experience™ is not the name of a function, but is Samsung's marketing terminology.
- \star S-Vision ${}^{\rm M}$ is the name of Samsung's ultrasound imaging technology.
- * S-Vue^ M is the name of Samsung's advanced transducer technology.
- * S-Detect™ for Breast and S-Detect™ for Thyroid are not available in Canada.
- * Strain value for ElastoScan+ $^{\rm TM}$ is not applicable in Canada and the United States.
- * Recommendations about whether results are benign or malignant are not applicable in the United States.

Empowering Professionals

Ultrasound System **RS85**



SAMSUNG MEDISON CO., LTD.

© 2017-2018 Samsung Medison All Rights Reserved. Samsung Medison reserves the right to modify the design, packaging, specifications, and features shown herein, without prior notice or obligation.

CT-RS85 V1.0_2-Eda-171201-EN

A New Healthcare Solution



Scan code or visit www.samsunghealthcare.com to learn more

SAMSUNG

A New and Outstanding






Improved Diagnostic Confidence

Samsung's image enhancing and artifact suppressing technologies and S-Vue™ transducers together provide clear, detailed imaging that you can count on to help improve diagnostic confidence and imaging continuity.

S-Vue[™] Transducers

S-Vue™ transducers provide more efficient piezoelectric properties, resulting in wider bandwidths that enable better penetration and higher quality resolution on even challenging patients.

 * Compared with the conventional Samsung transducers.
* The image is for illustrational purposes only and might differ from the actual performance of the device.



S-Vision[™] Imaging Engine

With the S-Vision[™] imaging engine built into RS85, the digital signals produce clear, detailed resolution and tissue uniformity for various types of applications in general imaging.

S-Harmonic[™]


This new harmonic technology improves image clarity, near to far. Reducing signal noise, S-Harmonic[™] provides more uniform ultrasound images. Combined with the S-Vue[™] transducers, S-Harmonic[™] takes RS85 image quality one step further.

HQ Vision™


HQ Vision™ is a new, advanced technology for visualizing anatomical structures. With improved image clarity, this feature helps make a reliable diagnosis quickly.

ClearVision


The noise reduction filter improves edge enhancement and creates sharper 2D images for optimal diagnostic performance. The integration of specialized Samsung technology results in a notable improvement in image quality. In addition, ClearVision provides applicationspecific optimization and advanced temporal resolution in live scan mode.




* The image above is for illustrational purposes only and might differ from the actual performance of the device.



Gallbladder



Finger



Pancreas



More Valuable Information

Expert tools offer new perspectives and provide additional information for confident decision making.

MV-Flow™

MV-Flow[™] offers a novel alternative to power Doppler for visualizing slow flow microvascularized structures. High frame rates and advanced filtering enable MV-Flow[™] to provide a detailed view of blood flow in relation to surrounding tissue or pathology with enhanced spatial resolution and temporal resolution.



Kidney

Thyroid

CEUS+

CEUS+ technology uses the unique properties of ultrasound contrast agents. When stimulated with low acoustic pressure, the oscillating microbubbles reflect both fundamental and harmonic frequency signals. In addition, Samsung's technologies, VesselMax[™] and FlowMax[™], provide a clear visualization of vessels and blood flow for a more informed and confident diagnosis.

S-Fusion[™]

S-Fusion™ enables simultaneous localization of a lesion using real-time ultrasound in conjunction with other volumetric imaging modalities. Samsung's Auto Registration helps quickly and precisely fuse the images, increasing efficiency and reducing procedure time. S-Fusion™ enables precise targeting during interventional and other advanced clinical procedures.

S-Fusion[™] for Prostate

S-Fusion™ for Prostate allows precise targeting during prostate biopsies. Based on 3D models created with MR data sets, S-Fusion™ for Prostate provides biopsy guidance to help safely navigate and target the prostate.




Increased Consistency

With its advanced intelligent solutions, including an extensive range of quantification functions, RS85 provides measurement consistency while reducing variability between users.

S-Shearwave Imaging™

S-Shearwave Imaging[™] allows for non-invasive assessment of the stiffness of tissue/lesions in the breast and liver, by providing an advanced level of diagnostic information. The color-coded elastogram, quantitative measurements (in kPa or m/s), dual or single display option, and user-selectable ROI (position and size) functions are especially useful for the accurate diagnosis of breast and liver diseases.



Arterial Analysis™

Arterial Analysis[™] detects functional changes of vessels, providing measurement values such as the stiffness, intima-media thickness and pulse wave velocity of the common carotid artery. Since the functional changes occur before morphological changes, this technology supports the early detection of cardiovascular disease.

S-3D Arterial Analysis™


S-3D Arterial Analysis™ simplifies volume measurement of arterial plaque, providing 3D vessel modeling. With Samsung's S-3D Arterial Analysis™, obtaining information on the arterial plaque volume is surprisingly fast and easy even on difficult patients. In addition, it allows you to track the morphological changes of the artery.

S-Detect[™] for Breast

S-Detect[™] for Breast helps standardize reporting and classification of suspicious breast lesions by incorporating BIRADS ® ATLAS* (Breast Imaging-Reporting and Data System, Atlas) into the tool. When the user selects a region of interest, S-Detect[™] for Breast automatically defines the lesion boundaries, provides lexicon classification options, and images export for an enhanced and streamlined workflow.

*Registered trademark of the American College of Radiology. All rights reserved.





Common carotid artery



S-Detect[™] for Thyroid

S-Detect[™] for Thyroid uses the advanced technology based on K-TIRADS, RUSS and ATA guideline* in detecting and classifying suspicious thyroid lesions semi-automatically. This technology helps you diagnose your patients with confidence and ease, providing accurate, consistent results and an automatic reporting feature.

* K-TIRADS: Korean-Thyroid Imaging Reporting and Data System RUSS: Russ' TIRADS ATA: American Thyroid Association




Enhanced Efficiency

The RS85 has been designed to streamline your workflow by enhancing efficiency through reducing keystrokes and by combining multiple actions into one.


Quick Preset

With one touch, the user can select the most common transducer and preset combinations. Quick Preset increases efficiency to make a full day of scanning simple and easy.



Touch Customization

Samsung has made a customizable touchscreen interface that allows the user to move frequently used functions to the first page, keeping the focus on the patient instead of the system.







6-way Control Panel

The RS85's 6-way adjustable control panel optimizes your work environment to reduce repetitive motions stress. When it's in off-mode, the control panel returns to the home position, allowing for easier and enhanced mobility.



13.3-inch Tilting Touch Screen —

Samsung's tilting touch screen can be adjusted to accommodate any user's viewing preferences within any scanning environment.

Central Lock

A single pedal controls a central lock mechanism to conveniently secure the console in place. This results in more efficient movements while the user is performing scanning procedures.



Streamlined Workflow



Gel Warmer

Samsung's two-level adjustable gel warmer keeps ultrasound gel at a comfortable temperature.



Maneuverable Wheel

4 swivel wheels allow easy steering, and a locking function.


SAMSUNG







Liver






E-Breast[™]








Patellar tendon in Panoramic




Breast mass




Internal carotid artery




Strain+



Common carotid artery



Auto IMT+



4 Chamber

Comprehensive Selection of Transducers



Linear array transducers









musculoskeletal

Application: small parts, vascular, musculoskeletal

Curved array transducers



Volume transducers



Phased array transducers



CW transducers



* Some of the transducers may not be available in some countries.





Application: small parts, vascular, musculoskeletal, abdomen

Application: musculoskeletal

L3-12A

LA3-16A

14 RS85

musculoskeletal





Application: cardiac, pediatric



Application: cardiac

TEE transducer



Application: cardiac